Carefully read the Instructions for Use. Observe all warnings and precautions. Failure to do so might result in complications.

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Five (5) MaRVis MR Guidewires in sterile bag per cardboard box.

Product lines for MaRVis MR Guidewires:

- MaRVis Amber Wire
- MaRVis Sunglow Wire
- MaRVis Tuscany Wire

Description

The MaRVis MR Guidewire is a medical device. It is a polymer-based, MR Safe guidewire that does not contain a metallic core. It is built from fiber-reinforced plastic rods as the base material. The guidewire diameter, length, guidewire type ("standard" or "stiff", if applicable) and flexible tip type (type A or type B, if applicable) are indicated on the product label. The flexible tip (distal part of brown section) is manually shapeable (only Type A, not Type B variants). The entire guidewire is covered by a PTFE shrink tube. The transparent PTFE surface of the MaRVis MR Guidewire may appear whitish in sections. This is a result of the manufacturing process and does not influence properties. The MaRVis MR Guidewire is visible under magnetic resonance imaging (MRI) over its entire length (continuous MR marker doping in the guidewire shaft) and has an additional MR tip marker (discrete larger artifact). For both markings, iron microparticles are used as MR markers. Only the tip (brown section) is visible under X-ray imaging.

Intended Use

The MaRVis MR Guidewires are intended for use in MRI-guided interventional medical procedures to reach the target region in the human body and/or to advance catheters or catheter-based instruments over the guidewire.

The MaRVis MR Guidewires must be used by a well-trained physician.

Scope of application for the MaRVis MR Guidewire product lines

A) MaRVis Amber Wire

Use up to one day in intraluminal interventions in natural body orifices in the urinary tract, the genital tract or the gastrointestinal tract and in non-vascular interventions in the trunk and the limbs.

The MaRVis Amber Wire is not intended for application in the peripheral or central circulatory system.

B) MaRVis Sunglow Wire

Use up to one day in endovascular interventions in the peripheral and central circulatory system. The MaRVis Sunglow Wire is not intended for application in the cerebrovasculature.

C) MaRVis Tuscany Wire

Use up to one day in endovascular interventions in the peripheral and central circulatory system.

Contraindication

Crossing of chronic total occlusions (CTO)

Warnings

- The MaRVis MR Guidewire does not contain a metal core. Therefore, it may present mechanical and handling properties differing from those of metal core guidewires.
- DO NOT CURVE THE SHAFT BELOW THE BENDING RADIUS AS INDICATED ON THE LABEL.
- Do not pull back when inserted into a metal needle, metal dilator, metal cannula, or metal introducer. Do not use with a sharp tool.
- Use caution when manipulating, advancing and/or withdrawing the MaRVis MR Guidewire through stents or other edged medical devices, or through tortuous or calcified blood vessels.
- Manipulation, advancement and/or withdrawal of the guidewire past sharp or beveled edges may result in destruction and/or separation in particular of the outer PTFE layer. This may lead to adverse clinical events.
- Do not apply local pressure, torsion, bending or transversal force. Do not apply local force to a particular point of the guidewire. This may cause damage to the guidewire. Do not affect the structural integrity of the guidewire.
- Replace any kinked, bent or otherwise damaged MaRVis MR Guidewire. It must not be used any longer because any damage to the guidewire may result in vessel damage and/or impaired performance characteristics. Fracture may lead to loss of torque control.
- When reintroducing the guidewire into the patient, it must be assured that no damage has occurred to the guidewire during preceding use. In case of detected or suspected damage do not use the guidewire anymore.

- When using the MaRVis MR Guidewire with a device that emits energy (heat, cold, laser, pressure, ultrasound, etc.) confirm that the MaRVis MR Guidewire will not be impacted by the energy.
- If using another interventional device comprising an active or resonant MR marker compatibility (e.g. interference between the markers) must be verified before use.
- The manually shapeable flexible tip (distal part of the brown section only!) may be shaped with fingers (with gloves!) but only cautiously. Do not use tools. Do not apply heat or steam. Attention: only Type A but not Type B variants are manually shapeable!
- A retrieving device such as e.g. a gripper, basket forceps, may only be used after the MaRVis MR Guidewire has been completely removed. Using retrieving devices while the MaRVis MR Guidewire is in the vessel or any other body tissue may cause the MaRVis MR Guidewire to break.
- Manipulate the MaRVis MR Guidewire slowly and carefully while confirming the behavior and the location of the tip in the anatomy under real-time image guidance. Excessive manipulation of the MaRVis MR Guidewire without confirmation by image guidance may result in vessel perforation or tissue damage.
- If any resistance is felt, or if the tip behavior and/or location seems improper, stop manipulating the MaRVis MR Guidewire and/or the catheter and determine the cause by appropriate imaging.
- If the MaRVis MR Guidewire is entrapped, e.g. in vasculature, do not torque, push or pull the MaRVis MR Guidewire. In appropriate cases, according to known medical retrieval techniques, placing a catheter over the guidewire for its retrieval may be useful.
- The introducer at the outer end of the dispenser tube is not intended to enter the body. It is only intended to facilitate insertion of the guidewire into e.g. hemostatic valves.

Precautions

- MaRVis MR Guidewires are delicate instruments and must be handled carefully. Failure to exercise proper caution during manipulation may result in bending, kinking, damage to the outer surface of the MaRVis MR Guidewire, damage to a catheter, or damage to the anatomy.
- The MaRVis MR Guidewire must be used by an experienced physician who is well trained in manipulation and observation of guidewires under real-time image guidance (MRI, fluoroscopy/angiography, as applicable).
- Advance the guidewire only a few centimeters at a time in the lumen of another medical device. If there is too much feed length between the hub (entrance) of the device and the location where the guidewire is grasped, the guidewire may kink. Reduce the feed length immediately, if the guidewire begins to bend.
- Procedures requiring the use of a guidewire shall not be attempted by physicians unfamiliar with the possible complications. **Possible complications** include, but are not limited to: Air embolism/thromboembolism, allergic reaction, amputation, arteriovenous fistula, burn, death, detachment of plaques, embolism, hematoma, hemorrhage, hemoglobinuria, infection or sepsis, irritation, inflammation, myocardial ischemia and/or infarction, pseudoaneurysm, stroke (CVA)/transient ischemic attack (TIA), thrombus, vessel occlusion, vessel perforation, dissection, trauma, damage, vessel spasm, guidewire entrapment, foreign body/guidewire fracture, ischemia.
- **Sterile.** The MaRVis MR Guidewire has been sterilized by ethylene oxide.
- Sterile in an unopened and undamaged sterile bag. Do not use if the sterile bag or the MaRVis MR Guidewire is broken or soiled. Do not use the MaRVis MR Guidewire if the guidewire is damaged.
- The MaRVis MR Guidewire should be used immediately after opening the package and be disposed of safely and properly after use, adhering to the local regulations for medical waste management.
- This medical product is designed for **single use only**. It must not be reprocessed, cleaned, disinfected, or sterilized. After reprocessing the product is no longer suitable for use and may constitute a biological risk for the patient and user. Due to the design of the product, efficacy of cleaning after use cannot be ensured, and contamination cannot be ruled out. Re-sterilization may affect the surface and impair material properties. In this case functionality of the product cannot be ensured.
- Consider using systematic heparinization.
- When using a drug or device with the MaRVis MR Guidewire, the operator should have full understanding of the properties/characteristics of the drug/device to avoid damage to the MaRVis MR Guidewire.
- Concurrent use of the MaRVis MR Guidewire with any other medical device must be verified before use. **The 0.035" MaRVis MR Guidewires are compatible with catheters of at least 5Fr size (inner diameter at least 1.00 mm).**
- Confirm the compatibility of the guidewire diameter with the other interventional devices before use.
- Flush the dispenser with heparin-saline solution before pulling the MaRVis MR Guidewire out of it. Prior to insertion of the guidewire into a catheter (or other device), flush the catheter with heparin-saline solution.

- After removal from the patient, and prior to reinsertion into the same patient during the same intervention, the MaRVis MR Guidewire shall be rinsed in heparin-saline solution.
- Use of alcohol, antiseptic solution, or other solvents is prohibited as this may result in damage to the MaRVis MR Guidewire.
- Never reinsert the MaRVis MR Guidewire into the dispenser as the edge of the dispenser may cause damage to the surface of the guidewire, or any other mechanical damage.
- Due to the slippery nature of the outer PTFE surface of the MaRVis MR Guidewire, the operator may encounter some difficulties in torquing the MaRVis MR Guidewire. An appropriate torque device is recommended for easier torquing of the guidewire. Do not use a metallic torque device as it may result in damage to the MaRVis MR Guidewire and may be unsafe in the MR environment.
- During insertion into a medical device damage to the surface of the MaRVis MR Guidewire may occur if this device has a too small inner diameter. If any resistance is felt during introduction into the lumen of a device it is recommended to stop using such device.
- Do not slip a tightened torque device or Y-connector over the guidewire or manipulate the MaRVis MR Guidewire through a tightened up rotating hemostasis valve as this may result in damage to the MaRVis MR Guidewire. Do not fix a torque device too strongly as it may damage the guidewire.
- Do not use after the use-by date has been reached.
- **Storage.** Always store in the white cardboard boxes and do not remove the sterile bags from cardboard box for storage. Keep away from sunlight, UV light, water. Keep dry. Do not expose to temperatures below 5° C or above 40° C.

In the MR environment

- The MaRVis MR Guidewire is MR Safe according to ASTM F2503-23.
- Visualization of the devices and of the anatomy under MRI-guidance is different from conventional fluoroscopy/angiography or CT imaging.
- The physician, as well as the interventional team, shall be well trained in the MR environment and skilled in MR imaging.
- The MaRVis MR Guidewire is visible under MRI over its entire length (continuous MR marker doping in the guidewire shaft). The additional MR tip marker (discrete larger artifact) at the very distal end of the guidewire indicates the position of the tip of the guidewire.
- The MR tip marker must be visible in the visualized MRI slice for correct localization of the distal tip of the MaRVis MR Guidewire throughout the procedure.
- Only appropriate MR sequence settings will provide satisfactory concurrent visualization of the body tissue and the MaRVis MR Guidewire as well as of potentially other devices used.

Directions for use

1. The MaRVis MR Guidewire is intended for use up to one day in contact with the human body.
2. Remove the dispenser containing the MaRVis MR Guidewire from the package without pulling the guidewire out of the dispenser.
3. Remove the stopper at the inner end of the dispenser tube, without pulling the guidewire out of the dispenser tube.
4. Fill the dispenser tube with sterile heparin-saline solution through the hub of the dispenser, using a syringe.
5. Remove the MaRVis MR Guidewire from the dispenser by carefully pulling the guidewire at the tip (brown section) until the guidewire shaft (yellow section) starts exiting the dispenser tube. Then grasp the shaft of the guidewire to completely pull it out of the dispenser.
6. Inspect the MaRVis MR Guidewire prior to use for presence of bends, kinks, surface irregularities, or damages. Do not use if any defects are detected.
7. Prior to use, rinse the catheter/interventional device with heparin-saline solution to ensure smooth movement of the MaRVis MR Guidewire within it.
8. Insert the flexible tip (brown section) of the MaRVis MR Guidewire ahead into the anatomy or the other interventional device used. Inserting the proximal end of the guidewire ahead may damage the anatomy.
9. Keep at least 5 cm of the MaRVis MR Guidewire protruding from the hub of the interventional device.
10. After use the MaRVis MR Guidewire can be disposed of with normal clinical waste.

Any serious incident that has occurred in relation to the MaRVis MR Guidewire should be reported to the manufacturer and/or its authorized representative/UK Responsible Person and/or the competent national authority of the country in which the user is established according to local medical device regulation.

Access to Eudamed databank: <https://ec.europa.eu/tools/eudamed/#/screen/home>

Liability and complaints

Liability. The MaRVis MR Guidewire is manufactured with the greatest possible care. MaRVis Interventional GmbH has no influence on how the product is used, the diagnostics or therapy of the patient, nor how the product is handled outside of the company. MaRVis cannot guarantee the effectiveness nor the absence of complications. Therefore, MaRVis assumes no liability for damages and costs. MaRVis will replace products that have a defect for which MaRVis is responsible. MaRVis cannot be held liable for consequential damages of any kind which was caused by non-observance of the aforementioned Instructions for Use (for example, re-sterilizing or reprocessing the product).

Complaints. In the event of a complaint, a reportable or potentially reportable incident, all products connected with the incident and made by MaRVis or any other brand must be preserved if possible for 3 months under the allowed storage conditions. If required, they must be handed over to the responsible authorities or, if necessary, to the manufacturer itself for further analysis.



Standards. Indications, contraindications, and application of the product must always comply with current medical standards. The guidelines and recommendations of the relevant specialist medical association must be observed.

MaRVis Interventional GmbH
Drosselstr. 3
82467 Garmisch-Partenkirchen
Germany
www.marvistech.com
info@marvistech.com

CE 0044

Symbols

	Year-Month	Country and date of manufacture
	Year-Month	Use-by date
		Do not curve the shaft below a bend radius of x mm!
		Store at temperatures between 5° and 40° C
		MR Safe
		Catalogue number
		Batch code
		Single sterile barrier system
		Single sterile barrier system with protective packaging outside
		Do not use if package is damaged
		Fragile, handle with care
		Keep away from sunlight
		Keep dry
		Sterilized using ethylene oxide
		For single use only, do not reuse
		Caution
		Medical device
		Manufacturer
		Consult Instructions for Use
QTY:		Quantity (number of MaRVis MR Guidewires in the pouch/cardboard box)
		CE Mark

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